


Automated Phantom Analysis



- Phantom templates with ROIs and built-in functions
- DICOM tag analysis to pre-process images
- Automatically inserts phantom and DICOM data into QC-Track worksheets
- More sensitive to subtle equipment performance changes than visual assessment
- Supports range of phantoms, including manufacturer, ACR/Gammex 464, and Catphan

Accurate: manual data entry is substantially reduced

Fast: phantom processing in seconds

Consistent: same analysis and scoring used for all imaging units

Powerful: multi-slice operations

Flexible: built-in editor makes it easy to make phantom-specific and site-specific changes

QC-Track: Complete Paperless QC Solution

With QC-Track, CT teams now have a system to help meet tracking, trending, and security requirements.

Communication


Emails keep the whole team aware of failed or overdue tests and warn if tests are approaching overdue status. Credential emails alert clinical staff to deficiencies and approaching renewal dates. Event and Document tracking centralize scheduling and emailing for physicist surveys, preventive maintenance, and service.

Accountability

QC Dashboards clarify expectations, highlight issues, and encourage proper QC workflows, including full-circle failure follow-up.

Security

On-prem packaging meets ever-stricter IT security requirements.



Interested?

Contact us to learn more about Atirix QC solutions.

877-273-1764 | Atirix.com | QC-Answers@atirix.com

Upgrade to an Enterprise QCIS™

Start with QC-Track for CT, and then enhance your investment and be Always Inspection Ready™ across modalities, including modules for Mammo/FDA EQUIP and Ultrasound HLD/Exam traceability.