


Module: Fluoroscopy QC

- Track key parameters, like kVp, mA, and dose index
- Improve fleet-wide consistency and QC sensitivity with optional phantom-based QC
- Store documentation logs with permanent retention
- Easily access physicist annual survey

Centrally Manage Fluoroscopy QC with QC-Track®




- Schedule workflows to gather performance data on all fluoro units, including interventional suites, cath labs, R/F rooms, C-arms, and mini C-arms
- Receive immediate email notifications of QC failures
- Monitor and manage QC with suite of reports
- Scan and attach documents such as Annual Survey reports and ongoing PM/Service reports
- Optional enhancements for a complete workflow, including the QC-Track Planner module to set up and track all fluoro QC events, including physicist surveys


Optional: Credential Tracking Module to Meet State, Certification, or Institutional Requirements

- Meet California Title 17 credential tracking requirement
- Track and report on the training of fluoro clinical professionals
- Automatically send email reminders and notices for expiring and overdue credentials
- Simplify compliance with institutional requirements by linking to training content directly from reminders

See our Credentials brochure for more information.



Research



Research demonstrated that signal-to-noise ratio is an effective and efficient quality measure for daily fluoroscopy QC. The QC-Track software with automated fluoroscopy phantom analysis allows both technologists and physicists to quickly obtain key QC values, making regular, image-based QC a realistic practice for busy IR labs.

Contact research@atirix.com for a copy.



Automated Phantom QC for Fluoroscopy



Leeds TOR 18FG

- Supports loop processing for QC reflective of clinical use
- Standardizes fleet-wide QC with consistent ROI placement and objective analysis
- Increases sensitivity to subtle equipment performance change
- Provides convenient analysis fast enough for daily QC
- Structures reports for ongoing workflow management, physicist oversight, and inspections
- Configurable for full range of phantoms and fluoro manufacturers

Recent Implementations Included:

- Siemens
- Philips
- OEC (GE)
- Shimadzu
- Orthoscan

Interested?

Contact us to learn more about Atirix QC solutions.

877-273-1764 | Atirix.com | QC-Answers@atirix.com

Upgrade to an Enterprise QCIS™

Start with QC-Track for fluoroscopy, and then enhance your investment and be Always Inspection Ready™ across modalities, including modules for Mammo/FDA EQUIP and Ultrasound HLD/Exam traceability.